



ECO2

Powerful and Affordable B&W system



CHISON
Value Beyond Imaging

User-friendly and Modern Design



12 inch rotatable LED monitor (0~30°)
with different viewing angles



One key to quickly save



Water-proof keyboard Cover



DICOM 3.0 (option), USB

ECO 2

Premium Design, Surprising Price



Ultra-compact



Clips for quickly recall
(CINE, Still images)



Carrying case (BG-100):
go anywhere you want

Advanced Imaging Technologies



Ergonomics

- Rotatable LED (0~30°): tailor to different viewing angles
- Ultra-compact & Ultra-light
- Easy to carry
- Keyboard: For easy and quick access
- Dedicated cart TR 9000
- Streamlined workflow

ECHO

- Advanced technologies:
 - i-Image (Intelligent image optimization)
 - THI (Tissue Harmonic Image)
 - SRA (Speckle Reduction Algorithm)
 - Compound imaging
- Full digital beam-former
- Trapezoid
- Automatic PW trace
- Chroma
- Full screen display mode
- AIO (Automatic image optimization)

Economic

- The latest platform & more advanced chip-set: brings you better image quality and features at less cost in a more compact box
- Support PC printer, reduce investment & day to day paper cost
- 8G memory card
- 500G harddisk (option)
- Vehicular Charging

ECO-friendly

- Less radiation
- Less power consumption
- Longer usage life: LED
- Longer battery > 2 hours

Support PC Printer



ECO2



Traditional Ultrasound

Connect PC Printer



PC Printer

A4 Paper



Benefits

- PC printer is less expensive than video printer
- A4 paper saves more cost than thermal paper for everyday use
- Detail report can be printed (with images, measurement results, patient information)

Connect Video Printer

Can ONLY Use Video Printer



Video Printer

Thermal Paper



- More expensive
- Can ONLY print images

Excellent PW resolution and automatic

PW trace measurement will greatly improve the diagnosis confidence by delivering both qualitative & quantitative results.

Great image quality



Median Sagittal Section,B Mode



Early Pregnancy,B Mode



Kidney,B Mode



Subscapularis Tendon,B Mode



Thyroid,B Mode



Gallbladder,B Mode



Aorta Artery,PW Mode



Thumb Artery,PW Mode



Carotid,PW Mode

Specifications

Applications	Abdomen, OB/GYN, Vascular, Small parts(breast, testicles, thyroid), Musculoskeletal Conventional & Superficial, Pediatrics, Basic Cardiac
Imaging technologies	i-Image, Compound imaging, SRA, THI
Display Modes	B, B/B, 4B, M, B/M, PW, CPA, DPD
Gray Scale	256 Levels
Scanning Depth	240mm Max(depend on the probe types)
Probe Frequency	2.5 MHz - 11.0 MHz (Depend on the probe)
Cine Loop	256 frames
Zoom	4 steps
Software Packages	Dedicated measurement package: General, OB/GYN, URO, Cardiac, Vascular, Small parts
Archives Management	Patient information management, images record, exam report
Power Supply	AC 100V-240V, 50Hz-60Hz
Standard Configuration	1 probe connector, 3 USB ports, VGA port, Video port, LAN port, Full Screen Display Mode, i-Image, Chroma, SRA, Compound imaging, THI, 8G memory card, PW Mode
Option	2 probe connectors, DICOM, Water-proof keyboard cover, Carrying Case (BG-100), Cart TR 9000, Support PC Printer, Video Printer, Built-in Battery



2.5MHz-5.0MHz Convex
C3-A



5.3MHz-10.0MHz Linear
L7M-A



5.3MHz-11.0MHz Linear
L7S-A



2.5MHz-5.0MHz Micro-Convex
MC3-A



5.0MHz-10.0MHz Trans-Rectal
R7-A



4.5MHz-8.0MHz Transvaginal
V6-A



4.5MHz-8.0MHz Micro-Convex
MC6-A

CHISON Medical Technologies Co., Ltd.

Sales & Service Contact Address: No.3, Changjiang South Road, Xinwu District, Wuxi, Jiangsu, China 214028

TEL : 0086-510-85310937 **FAX :** 0086-510-85310726 **EMAIL:** export@chison.com.cn

We reserve the right to make changes to this catalogue without prior notice. All rights reserved.
Please contact our local dealer for the latest information.

Profimédica
INNOVANDO EL FUTURO



www.profimedicave.com



profimédica



04129090462